

Spectrum of efficacy of our product Endo Star

Endo Star is a quaternary ammonium compound- and phenol-free disinfectant concentrate for the disinfection of flexible and rigid endoscopes, endoscope accessories and anaesthetic accessories.

It is a medical product of the class IIa and is CE 0481 marked.

It is registered with the DGHM list of disinfectants and has the corresponding certificate dated of 31 December 2003 with the following values:

1 %	-	60 minutes
2 %	-	30 minutes
3 %	-	15 minutes

(a new and updated VAH/DGHM certificate is under way and supposed to be available in a couple of weeks)


Endo Star has the following spectrum of efficacy:

- **bactericidal, fungicidal, tuberculocidal, virus inactivating**

Successfully tested on the following disease causing agents:

- **HBV, HIV, BVDV-HCV, Vaccinia, Staphylococcus aureus, Escherichia coli, Proteus mirabilis, Pseudomonas aeruginosa, Candida albicans, Enterococcus hirae, Mycobacterium terrae, Aspergillus niger.**

Kempen, the 24th of April 2007


LABORATORIUM
Dr. H. D. Deppe
Hooghe Weg 35
47906 Kempen
Tel. (0 21 52) 55 65 -0
Fax (0 21 52) 5 08 49
Titus Deppe
Hygienic Products Manager