

Top Off wipes

alcoholic, aldehyde free

Medical product of the class IIa, CE 0481

Article number 90200 (150 wipes box), 90201 (150 refill-package)

Type of product:

Hygienic ready-to-use wipes (soaked with SPRAY OFF N) for the cleaning and disinfection of all alcohol resistant surfaces and equipments, also in the proximity of patients. They are very economical due to the use of re-fill-boxes. The wipes are lint-free, not sticky, do not dry up and leave no visible residue. They are perfect for small surfaces and small equipments.

The unfolded size of each wipe is 13,5 x 20 cm.

Note: Be careful with materials which are sensitive to alcohols and with alcohol-dissoluble paint coats and varnishes.

Storage conditions: do not allow temperatures below 0 °C or above 30 °C; do the product not expose to direct sunlight; keep the store dry



Surfaces

Surfaces cleansing and disinfection



Ingredients and spectrum of effect:

see SPRAY OFF N.

Main advantages:

+ highly effective; + quickness; + pleasant smell; + gentle on instruments; + low alcohol contents; + aldehyde free; + good skin tolerance with short contact; lint free

Available sizes:

150 wipes box
150 wipes re-fill-package

Shelf life:

3 years

Top Off XXL wipes Top Off XXL Alpha wipes

alcohol based (Spray Off N)

alcohol free (Alpha Guard Sensitive GF)

Medical product of the class IIa, CE 0481

Article number

90202 Alpha Guard GF Sensitive Wipes (1 bucket + pouch x 70 wipes)

90209 Alpha Guard GF Sensitive Wipes (6 refill pouches x 70 wipes)

90203 Top Off XXL (1 bucket + pouch x 70 wipes)

90204 Top Off XXL (6 refill pouches x 70 wipes)

Type of product:

- hygienic ready-to-use wipes for the cleaning and disinfection of surfaces and equipments
- for alcohol sensitive surfaces use the Alpha Guard GF Sensitive Wipes (1 min.)
- for alcohol resistant surfaces use the Top Off XXL wipes (30 sec.)
- size of the wipes: 26 x 29cm



Ingredients and spectrum of effect:

Top Off XXL: see Spray Off N

Alpha Guard GF Sensitive:

see Alpha Guard GF

Main advantages:

+ highly effective; + quickness (Top Off XXL: 30 sec; Alpha Guard GF Sensitive: 1 min.) + pleasant smell, + gentle on instruments; + aldehyde free (Top Off XXL); + alcohol and aldehyde free (Alpha Guard GF Sensitive); + good skin tolerance with short contact time; + lint free

Available sizes:

1 bucket + 1 pouch x 70 wipes
6 refill pouches x 70 wipes

Shelf life:

3 years