

**Dr. Deppe** 德國德布博士 德國原裝進口

## Thixo - Gel -Spray (ready-to-use) 超音波診療噴霧凝膠

Ultra sonic contact spray for the sonography (超音波診療專用)  
Aerosol(氣溶膠)

德國原裝進口，不需泡製，直接使用。

### Type of product:

Thixo Gel Spray is a universally applicable contact gel for the medical sonography. Due to the unique thixotropic properties THIXO GEL SPRAY can be easily sprayed and forms directly on the skin a viscous gel consistency.

超音波診療專用噴霧凝膠，使用方便，直接噴灑在受診部位，呈現液體凝膠，利於超波診療。

THIXO GEL SPRAY is pleasant for the application and does not glue and can easily be removed.

使用舒適簡便，不會膠黏皮膚感覺，擦拭容易。

- THIXO GEL SPRAY is based on mineral gel forming agents  
主要以礦物膠凝體成份形成。
- free of oils, polyacrylates and perfumes 非油性物質，芳香天然樹脂。
- highly conductive 高傳導性。
- precise contrast reproduction。高精確度光滑和柔順的圖像再生。
- with the skin friendly extra additive ALOE VERA。  
填加蘆薈增加與皮膚的親和力。
- quick and large area covering of the to be examined body parts。  
受檢身體部位能快速及大面積覆蓋。
- no necessity to spread the gel afterwards 不需要再使用其他凝膠。

### Application: 使用方式

Shake container well before usage and spray then the gel onto the to be examined body parts

使用前搖動(搖一搖)瓶身，噴口直接噴灑受檢部位。

**Storage conditions:** do not allow temperatures below 0°C or above 30°C;  
do the product not expose to direct sunlight; keep the store dry.

儲存方式：溫度 0~30°C，避免日光照射太久，維持產品乾爽環境。

skin



包裝：  
90115(250ml/瓶含噴頭，12瓶/case)

### Composition : 成份

Water, silicates, Aloe Vera extract, citric acid; Preservatives: Ethylparaben, 2-Bromo-2-Nitropropane-1,3-Diol, Cetrimonium Bromide. Propellant: Dimethylether

### 適用場所：

超波診療室

**Shelf life(儲存時間): 3years (出廠起3年)**

 **廣林貿易有限公司**  
KWANG LIN SCIENTIFIC TRADING CO., LTD.

高雄市80251民權一路34號14樓(超群大樓)  
14FL No. 34 Min Chuan 1st Rd. 80251 Kaohsiung TAIWAN

Tel: +886.7.3333485 Fax: +886.7.3366016

Kwanglin@ms26.hinet.net