

# Thixo - Gel - Spray

Ultra sonic contact spray for the sonography

Medical product of the class I article number 90115 (250ml bottle)

## Type of product:

Thixo Gel Spray is a universally applicable contact gel for the medical sonography.

Due to the unique thixotropic properties THIXO GEL SPRAY can be easily sprayed and forms directly on the skin a viscous gel consistency. THIXO GEL SPRAY is pleasant for the application and does not glue and can easily be removed.

- THIXO GEL SPRAY is based on mineral gel forming agents
- free of oils, polyacrylates and perfumes
- highly conductive
- precise contrast reproduction
- with the skin friendly extra additive ALOE VERA
- quick and large area covering of the to be examined body parts
- no or almost no necessity to clean up with paper

## Application:

Shake the container and spray the gel onto spots of the to be examined body parts and spread from there.

## Storage conditions:

do not allow temperatures below 0 °C or above 30 °C ; do the product not expose to direct sunlight; keep the store dry.

## Skin



## Composition:


Aqua, Glycerin, Propylene Glycol, Benzyl Alcohol, Methylchloroisothiazolinone, Methylisothiazolinone, Lithium Magnesium Sodium Silicate, Citric Acid, Aloe Vera Barbadosensis.

## Available Sizes:

250ml pump spray bottle

## Shelf life:

3 years

  
**THIXO-GEL ULTRASOUND SPRAY**  
Ultra sonic contact spray for the sonography

**Type of product:** Thixo Gel Ultrasound Spray is a universally applicable contact gel for the medical sonography. Due to the unique thixotropic properties Thixo Gel Ultrasound Spray can be easily sprayed and forms directly on the skin a viscous gel consistency.


Thixo Gel Ultrasound Spray is pleasant for the application, leaves no sticky residue and is very easily and quickly removed.


- Thixo Gel Ultrasound Spray is based on mineral gel forming agents
- free of oils, polyacrylates and perfumes
- highly conductive
- precise contrast reproduction
- with the skin friendly extra additive ALOE VERA
- quick and large area covering of the to be examined body parts
- no necessity to spread the gel afterwards
- dermatologically tested – harmless to skin and ultrasonic probe


**Application:**  
Spray the gel on a certain spot and distribute from there.

**Inhaltsstoffe:**  
Aqua, Glycerin, Propylene Glycol, Benzyl Alcohol, Methylchloroisothiazolinone, Methylisothiazolinone, Lithium Magnesium Sodium Silicate, Citric Acid, Aloe Vera Barbadosensis.

**Attention:** Application only externally!







Laboratorium  
Dr. rer. nat. Deppe  
D-47906 Kempen  
Germany  
www.dr-deppe.de

Content: **250 ml**

120510001