



TOTAL AEROBIC BACTERIA for sanitary monitoring
 總菌數培養片 (1.0mL培養後出現紅色菌落表示陽性反應)

Incubation at 35 ~ 37 °C for 24 ~ 48 hour

Catalogue number: BAC-096




Total Yeast & MOLD for sanitary monitoring

酵母菌與霉菌培養片 (1.0mL培養後出現藍色菌落表示陽性反應)

Incubation at 28 ~ 30 °C for 48 ~ 72 hour

Catalogue number: YFC-096



Total COLIFORM for sanitary monitoring

大腸桿菌群培養片 (1.0mL培養後出現紅色菌落黃色底部背景表示陽性反應)

Incubation at 35 ~ 37 °C for 15 ~ 24 hour

Catalogue number: CCL-096



Total COLIFORM & E. Coli for sanitary monitoring

大腸菌群與大腸桿菌培養片

(1.0mL培養後出現紅色菌落表示大腸桿菌陽性反應)

(培養後出現藍色菌落表示大腸桿菌陽性反應)

Incubation at 35 ~ 37 °C for 15 ~ 24 hour

Catalogue number: ECL- 096



E. COLI O157:H7 & E. COLI for contamination detection

大腸桿菌O157:H7快速培養鑑定片

(0.5mL培養後出現紅色菌落表示 E.Colio157 陽性反應)

(培養後出現藍色菌落表示 E.Colio 陽性反應)

Catalogue number: EOH-096



SALMONELLA for contamination detection

沙門氏菌培養片

(1.0mL培養後出現深褐色菌落表示陽性反應)

Incubation at 35 ~ 37 °C for 24 ~ 48 hour

Catalogue number: SAL-096